

Medela Clario suction pump

Pump Specifications	
Item numbers	AC-version
	014.0010 100–240V/50–60Hz Euro plug 2xA
	014.0011 100–240V/50–60Hz Euro plug A
	014.0014 100–240V/50–60Hz Euro plug B
	014.0040 100–240V/50–60Hz UK Plug A
	014.0050 100–240V/50–60Hz AUS/NZ plug A
	014.0020 100–240V/50–60Hz US plug A
	AC/DC version
	014.0110 100–240V/50–60Hz Euro plug 2xA
	014.0111 100–240V/50–60Hz Euro plug A
	014.0114 100–240V/50–60Hz Euro plug B
	014.0140 100–240V/50–60Hz UK Plug A
	014.0150 100–240V/50–60Hz AUS/NZ plug A
	014.0120 100–240V/50–60Hz US plug A
	Deluxe version (AC/DC + SoftCase)
	014.0210 100–240V/50–60Hz Euro plug 2xA
	014.0211 100–240V/50–60Hz Euro plug A
	014.0214 100–240V/50–60Hz Euro plug B
	014.0240 100–240V/50–60Hz UK Plug A
	014.0250 100–240V/50–60Hz AUS/NZ plug A
014.0230 100–240V/50–60Hz US plug A	
Instructions for use	177.1405 EN, DE, FR, IT, NL, SV, FI, NO, ES, PT (A) 177.1406 ET, HU, CS, LT, LV, PL, SL, SK, DA, EL (B)
Service manual	200.4001 EN, DE, FR
Max. vacuum [kPa/mmHg] Measured at 0 m, atmospheric pressure: 1013.25 hPa Please note: vacuum levels may vary depending on location (meters above sea level, atmospheric pressure and temperature).	Tolerance: min. –18/–135 ± 15% med. –36/–270 ± 15% max. –80/–600 – 15%
Flow [l/min]	15 l/min. ± 10%
Dimensions [mm]	223 x 255 x 95
Weight [kg]	AC 1.7 AC/DC 2.0
Drive unit	Membrane
Power supply	100–240V, 50–60 Hz
Power consumption [V]	12V, 3.4 A
Noise level [dB(A)]; measured in Medela acoustic room at 1 m distance	64.3
Protection against electric shock	Type BF
Medical Device Classification – According to 93/42/EEC	Class IIa
Protection against electric shock	Class II Equipment
Ingress Protection	IP21
Battery	NiMH
Operation from battery (at full charge)	>50 minutes
Charging time of battery	5 hours

Overflow protection	Threefold hygienic design: – Floater in canister prevents overflow – SafetyChamber – Membrane system No overflow can occur into the pump
Warranty [years]	2; Battery 6 months
Operating temperature, relative humidity, pressure	5 to 40 °C, 15–93 %, 70–106 kPa
Transport/storage temperature, relative humidity, pressure	–20 to +50 °C, 15–93 %, 70–106 kPa
Made in	Switzerland

Accessories	Item number
Canister set	074.0001
SafetyChamber	074.0002
Disposable Fingertip tubing	074.0003
Car connection cable	074.0017
SoftCase carrying bag	074.0007
Vacuum gauge	074.0006

CAUTION: Federal US law restricts this device to sale by or on the order of a physician.

Medical Vacuum Technology
for Healthcare Professionals

Please contact us or your local Medela
representative for details.

 Medela AG
Lättichstrasse 4b
6341 Baar, Switzerland
www.medela.com

Local contact: